



Abstract

- In 2001, Burnett, Kim & Lane-Smith [1] introduced the technique for measuring radon in water using a Durrige RAD7 continuous radon monitor [2] paired with a continuous water-air equilibrator, commercialised as the RAD Aqua. The RAD H2O, for discrete measurements, followed several years later [3].
- In the two decades since, these techniques have gained popularity for their ease of use and fast results, and radon measurements with the RAD Aqua and RAD H2O are now routinely used by aquatic scientists to quantitatively study groundwater movement in coastal, lacustrine and riverine systems.
- To date, several thousand scientific works have been published using these devices and techniques pioneered by Durrige.
- Here, we present the RAD7's successor: the RAD8 continuous radon monitor. Redesigned from the ground up with aquatic science applications in mind, the RAD8 represents a technological leap forward in the researcher's toolkit.

70% Higher Sensitivity

- RAD8's patent-pending measurement chamber geometry has been optimised using several approaches, including full 3D electrostatic field analyses on a per-volume basis and Monte Carlo simulations, with models validated through real-world testing.
- Commissioning tests and calibration data show that RAD8 achieves a sensitivity 70% greater than that of the RAD7, despite having a smaller measurement chamber.

Fully Ruggedised and Waterproof (Even with the Lid Open!)

- RAD8 is IP67 waterproof and dust proof, and designed to float.
- Can be operated in wet or dusty environments, or even at sea, without fear of damage from incidental water contact.
- Dual controls: touchscreen and physical membrane buttons.




Benchmarking Test: MnO2 Fibre with radium-226 & thorium-232

- Henrietta Dulai – University of Hawai'i.
- Using the method first described in [4].
- RAD7: 104.7 +/- 8.6 Bq/m³**
- RAD8: 102.5 +/- 15.4 Bq/m³**
- RAD7: 37.5 +/- 1.8 Bq/m³ (thoron)**
- RAD8: 37.5 +/- 1.7 Bq/m³ (thoron)**

Benchmarking Test: RAD H2O

- Measuring a radium-free radon-in-water source [5] in a 40 ml RAD H2O sample vial.
- Reference: RAD7 running in built-in WAT40 mode.
- RAD8 uses calculated air-water conversion factors, rather than the empirically determined ones that were used by the RAD7's WAT40 & WAT250 modes [6]. This brings RAD H2O in-line with Capture's Big Bottle analysis tools.

- RAD7: 83.46 +/- 8.14 Bq/L**
- RAD8: 86.4 +/- 2.33 Bq/L**





Wireless Connectivity & Enhanced Data Storage

- Built-in WiFi connectivity, as well as one COM port (data download and remote control) and two USB ports for accessories.
- Onboard storage for orders of magnitude more data than the RAD7, including full-resolution spectral information in 5-minute 'time slices'.

Benchmarking Test: RAD Aqua

- Henrietta Dulai – University of Hawai'i
- RAD Aqua water-air exchange setup of [1], modified to incorporate RAD8 and RAD7 in series, sampling air in a closed loop with the RAD Aqua.
- Response to low-high-low step-changes in radon concentration are shown below.
- Both response time and radon concentration in very good agreement with RAD7 reference.

- RAD7: 24.9 +/- 3.3 Bq/L**
- RAD8: 24.6 +/- 3.0 Bq/L**




On-Device Data Visualisation

- View radon / thoron time series data on-screen as it comes in, to get a picture of the trend over time at a glance.
- Spectrum display empowers users to make their own data quality checks in real time, allowing them to identify problems and fix them in the field.
- Avoid the need for repeat measurements!
- RAD8 stores both reduced data and full-spectrum information in 5-minute 'time slices', which can be combined and split in Capture after the event – never worry about setting the wrong Cycle length again!


Measure up to 2,500,000 Bq/m³

- RAD7 shaping time: 100 μ s. Digitiser frequency: 40 kHz.
- RAD8 shaping time: 4 μ s. Digitiser frequency: 30 MHz.
- This speed increase enables the RAD8 to measure extremely high radon / thoron concentrations before pulse pileup and deadtime begin to impact linearity.



Real-time readout


- RAD8 displays continuously updating radon and thoron activity concentrations, with their associated uncertainties.
- Users thus have access to the most up-to-date information, and an early warning of changing radon / thoron activity concentrations.



Benchmarking Test: Custom Big Bottle Setup

- Joseph Tamborski – Old Dominion University.
- Water samples collected from the Elizabeth River in 8L jerricans.
- Custom setup akin to RAD H2O & Big Bottle System.
- Water sample never comes into contact with air, which eliminates a major potential source of sampling error.

- RAD7: 10.67 +/- 1.15 dpm/L**
- RAD8: 10.95 +/- 1.21 dpm/L**



Capture Pro 8 with Capture Cloud

- Capture Pro 8 software for Windows and macOS offers an expanding suite of data graphing and analysis tools.
- Capture Cloud provides automated uploading and retrieval of RAD8 data, and the sharing of data with colleagues around the world.
- Capture Cloud Data Browser offers powerful tagging and search functionality, making it easy to organize and locate RAD8 data from any project, device, date range, or location.

References

[1] Burnett, W. C., Kim, G., & Lane-Smith, D. (2001). A continuous monitor for assessment of ²²²Rn in the coastal ocean. *Journal of Radioanalytical and Nuclear Chemistry*, 249(1), 167–172. <https://doi.org/10.1023/A:1013217821419>

[2] Durrige Company Inc. RAD7 Radon Detector <https://durrige.com/products/rad7-radon-detector/>.

[3] Lee, J. M., & Kim, G. (2006). A simple and rapid method for analyzing radon in coastal and ground waters using a radon-in-air monitor. *Journal of Environmental Radioactivity*, 89(3), 219–228. <https://doi.org/10.1016/j.jenvrad.2006.05.006>

[4] Kim, G., Burnett, W. C., Dulajova, H., Swarzenski, P. W., & Moore, W. S. (2001). Measurement of ²²⁴Ra and ²²⁶Ra Activities in Natural Waters Using a Radon-in-Air Monitor. *Environmental Science & Technology*, 35(23), 4680–4683.

[5] Kitto, M. E., Menia, T. A., Bari, A., Fielman, E. M., & Haines, D. K. (2010). Development and intercomparison of a reusable radon-in-water standard. *Radiation Measurements*, 45(2), 231–233. <https://doi.org/10.1016/j.radmeas.2010.02.002>

[6] Lee, K., & Burnett, W., Determination of air-loop volume and radon partition coefficient for measuring radon in water sample, *J. Radioanal. Nucl. Chem.* 298 2 (2013) 1359–1365.

Acknowledgements

- Thanks to Henrietta Dulai (University of Hawai'i) for benchmarking data, beta testing and photos from testing in Rapa Nui, including the background image.
- Thanks to Joseph Tamborski (Old Dominion University) for benchmarking data, beta testing, and photos from the Elizabeth River RAD7-RAD8 comparison.
- Thanks to Chris Smith (USGS) for detailed information about the protocol used to collect the Elizabeth River data.